



# Application Driven and Easy to Operate Sample Introduction Solutions



**Dr. Maja Budanovic**

**ICP Product Manager**

**Glass Expansion GmbH**

**Contact: [mbudanovic@geicp.com](mailto:mbudanovic@geicp.com)**

# Who is Glass Expansion?

GE have been specializing in sample introduction components from the probe to the cones — for ICPs since the early 1980s

- **Many ICP vendors package GE parts as part of the standard configuration.**
- **Manufacturers Supported:** Thermo®, Agilent®, PerkinElmer®, Shimadzu®, Analytik Jena®, Spectro™, Others

**Products offered:**



# Is Your Sample Introduction Optimized for Your Application?

Your application works best with a tailored solution — not a default setting.

## Optimization depends on your priorities:

- Maximum sensitivity
- Improved precision and reproducibility
- Robustness for high-matrix samples
- Minimal carryover
- Faster washout for high throughput
- Low-volume or low-flow sample compatibility
- Compatibility with challenging acids or solvents

**Note:** 99% of analytical problems occur within the sample introduction configuration.


## Application Dedicated SIS solutions:

[www.geicp.com](http://www.geicp.com)

Click on: [View Recommended Products for your Application](#)

### Applications

- [Animal feed](#)
- [Brines and salts](#)
- [Chemicals and fertilizers](#)
- [Clinical and forensic materials](#)
- [Drinking, ground and surface water](#)
- [Food and drink](#)
- [Geological with HF](#)
- [Geological without HF](#)
- [Isotopic Analysis of Minerals](#)
- [Metals](#)
- [Petrochemicals](#)
- [Plants](#)
- [Semiconductors](#)
- [Soil and sediment with HF](#)
- [Soil and sediment without HF](#)
- [Waste water and sludge](#)
- [Wear Metals in oil](#)






# Nebulizer Selection

Selecting the right nebulizer requires careful consideration of various factors:

Nebulizer		TDS (%)	Particulates (µm)	HF	Precision	Purity	Material
SeaSpray™		20	75	No	High	Good	Glass
MicroMist™		15	40*	No	High	Good	Glass
Conikal™		5	75	No	High	Good	Glass
Slurry™		1	150	No	High	Good	Glass
Quartz SeaSpray™		20	75	No	High	Excellent	Quartz
OpalMist™		15	75*	Yes	High	Excellent	PFA
DuraMist™		30	75*	Yes	High	Good	PEEK
VeeSpray™		30	300	Yes	Moderate	Good	Ceramic

\* Varies with nebulizer uptake





## Important Nebulizer Operating Parameters

Example: GE P/N A13-[1](#)-[UM04](#)

- Optimum nebulizer gas flow = 1.0 L/min (40 psi)
- Sample uptake rate ≤ 0.4 mL/min


# Optimal Nebulizer

## 1. Solutions Without HF




# Optimal Nebulizer

## 2. Solutions Containing HF




## 3. Organic Solutions



[Click Here - Care of Nebulizers](#)



# Precision & Sensitivity: Design Considerations

Quality of Aerosol  $\propto$  Quality of Results



Smaller Droplets Require Less Energy = Efficient Ionization

# Spray Chambers: Selection






[Click Here - Spray Chamber Solutions for Organic Applications](#)



# Spray Chambers: Helix CT Interface

Helix CT: Constant Torque = Reproducible day-to-day ICP Performance



[Download the Helix CT ICP Spray Chamber Application Note](#)

[Click Here - Care of Spray Chambers](#)

# Tracey™ BC Spray Chamber: Design Considerations



## Challenges in Routine ICP Work

- Frequent O-ring wear and replacement
- Long washout times and memory effects
- Poor wetting or carryover with HF or harsh matrices

## How the BC Design Helps

- **No O-Rings** → Less maintenance, faster washout
- **Low-Volume 30 mL Cyclonic** → Shorter stabilization, better throughput
- **Broad Compatibility** → Fits E-Torch, D-Torch, SDT/FDT
- **Cost-Effective** → Practical choice for routine analysis
- **Improved reproducibility** in maintaining tighter overall size specifications

## Washout Profiles for 1 ppm Hg



**Tracey™ Ball Joint Connection  
(BC) Spray Chambers**

# Exploring the Tracey™ BC PEEK Spray Chamber for Your Workflow

- **PEEK Construction**→ Good chemical resistance (up to 5% HF)
- **Superior Wetting:** PEEK material maintains excellent wetting properties with routine laboratory cleaning.
- **No Internal Surface Treatment:** Unlike TFE or PFA, this spray chamber requires no internal surface treatment.

## Comparison of Tracey BC PEEK to the PTFE Tracey

Below are the average intensity and RSD results from 41 optimization checks using the PEEK and PTFE spray chambers.

- The Optimization Solution contains 2 ppm Pb, As, and Mn in 1% HNO<sub>3</sub>.



	% Increase in intensity	%RSD
Pb	74%	0.65
As	90%	0.69
Mn	68%	0.68



*\*Comparison conducted by Specialty Chemicals Manufacturer - USA*

# Aggressive Sample Matrices: Torch Selection

**Examples:** Soils, wastewater, brines, high-acid digests, organics, lithium fusions

**Challenges:**

- High salt deposits and plasma temperatures shorten quartz torch life
- Frequent torch replacement increases cost of ownership
- O-ring failures or gas leaks can destabilize plasma or prevent ignition


**Mitigation:**

- Use a Demountable Torch → replace only the outer tube, lowering cost of ownership
  - *Narrow bore quartz: 1.0mm or less:* volatile organics
  - *Large bore quartz: 2.0mm or greater:* High TDS
  - *Ceramic (alumina):* HF-containing samples
  - *Platinum/Sapphire Injectors:* Inert applications
- Ferrule-based design: Secure injector seating, fewer leak points




**Upgrade:** Optional ceramic outer tube → resists devitrification, lasts longer, and maintains plasma stability

**Benefits:** Ideal for high-TDS, salty, or organic samples; hotter, more robust plasma improves sensitivity



P/N 30-808-4388  
E-Torch for Thermo® PRO Duo

Comparison of Quartz Tube set to Ceramic Outer tube set		
Element	% Increase in Sensitivity	%RSD
Zn (213) λ	17%	0.36
Ni (231) λ	19%	0.57
Mn (257) λ	14%	0.52



Ceramic Outer Tube Set  
P/N 31-808-4502

# Samples That Contain Particulates (1/2):


## Why address particulates issue?

- Particulates can clog fine sample lines or nebulizer channels, disrupting analysis
- Preventing blockages reduces downtime and maintains performance

## How to solve it?

**1a. Guardian** In-Line Particle Filter (P/N [70-803-1108](#)) between probe and nebulizer

- **120µm filter** with seals for both 1.6mm and 1.3mm OD tubing
- **Clog-resistant PEEK** design — easily cleaned by back-flushing or ultrasonics



**In-Line particle filter:** *“By the way, the particle filters that we have purchased are working out very well with our soil sample analyses on our ICP-OES units, have saved a lot of headaches with blocked nebulisers!”*

**Soil & Plant Laboratory - Australia**


## 1b. Guardian Sample Probe




# Samples That Contain Particulates (2/2):

## Guardian Sample Probe:

- Proprietary mechanical finish for superior wetting characteristics
- **Drip-resistant** to minimize cross-contamination, especially with oils.
- **Unique inbuilt particle filtering** prevents blockages in your nebulizer and capillary tubing.
- **Completely inert construction** (Ceramic, PEEK, and PTFE) for strong acid/solvent resistance.
- **PEEK sheath** designed to ensure precise alignment within the middle of the vial every time.
- **Interchangeable UniFit sample lines** (3000mm in length) to accommodate various IDs (e.g. 0.3, 0.50, 0.75 & 1.0mm)



Guardian Probe Assembly  
for ASX-200, 500, 800 Series  
P/N 70-803-1787



# Samples With High TDS: Argon Humidifier

## Why address high TDS issues?

- Salt deposits form at the nebulizer and injector tip
- Leads to analytical drift or even plasma extinguishing

## How to solve it?

### 1. Elegra™ Argon Humidifier to prevent salt build-up

- Adds moisture to argon gas → prevents salt build-up and reduces maintenance
- Maintains stable plasma and consistent results during high TDS analysis
- **Elegra:** No power, heat, or pressurization required – compact, inert, and easy to integrate
- Flexible configurations: single- or dual-channel versions with custom gas fittings
- Complements a high-TDS sample intro setup: **SeaSpray™ / DuraMist™ nebulizer, Twister™ spray chamber, and wide-bore injector**

(In reference to the Elegra: “Talking with my operators that are here today neither of them has changed a nebulizer since we put it on... We had been replacing nebulizers after about a week and a half... I will be ordering 2 more.

Contract Laboratory - USA



SeaSpray™





DuraMist™


# Samples With High TDS: Argon Humidifier

## Performance

- Includes easy-use bypass switch to disable humidification without disconnecting lines
- Tested to deliver over **4 hours** of stable internal standard signal with high-salt samples, outperforming other humidifiers
- Superior performance: up to 60% more effective relative humidity than competing models



## Elegra Application Note



# General Guidelines On Cone Materials

## Nickel Cones:

- **Balanced cost & performance;** standard for many applications
- Good **thermal & chemical resistance;** less prone to corrosion and deposition
- Runs **hotter** than copper, stays cleaner longer, more stable signals
- Suitable **for routine aqueous samples** (<5% acid, non-HF, non-organic)



## Nickel-Plated Cones:

- Ideal for samples with >5% acid concentration
- Nickel plating boosts chemical resistance while **retaining copper's efficient heat transfer**
- **Helps prevent overheating** and rapid orifice degradation, preserving sensitivity and stability

## ICP-MS Cone Resource Guide



## Platinum Cones:

- Most durable, longest-lasting, but highest cost
- Excellent chemical resistance - ideal for **high-matrix, high-acid**, or organic **solvent samples**
- **Least efficient** heat transfer → runs **hotter**, but stays **cleaner longer**
- Can be **refurbished 2-3 times** and **recycled** for reclaim value towards future purchases



# Summary: Steps to Overcome SIS Challenges

## 1. Improve Data Quality:

- Select appropriate nebulizer, spray chamber, torch/injector, and cones
- Tailor components to sample type for accuracy, precision & sensitivity

## 2. Maximize Sample Throughput:



- Address carryover and washout issues to improve efficiency

## 3. Enhance Performance with Accessories

- Use tools like Elegra, Trident CT, or Guardian In-Line Filter to improve stability

## 4. Ensure Longevity & Consistency

- Implement proper care and cleaning routines
- Reduce downtime through preventive maintenance



[Click here to view our Product Care Page](#)





# Thank You



## Glass Expansion - Europe Weilburg, Germany

### Europe

Friedenbachstrasse 9  
35781 Weilburg  
Germany

**Phone:** +49 6471 3778517

**Email:** gegmbh@geicp.com

**[www.geicp.com](http://www.geicp.com)**