

From Seawater to Shellfish: Microplastics... Find out what's slowly krilling you, and the best way to stay happy as a clam!



Ryan Brennan
President
Glass Expansion, Inc.
Rbrennan@geicp.com



Aaron Hineman
Inorganic Product Line Leader,
Americas
PerkinElmer
Aaron.Hineman@perkinelmer.com



GLASS EXPANSION
Quality By Design



Glass Expansion Offices

Head Office

6 Central Boulevard
Port Melbourne VIC 3207
Australia

Phone: +61 3 9320 1111

Email: enquiries@geicp.com

Americas

31 Jonathan Bourne Drive,
Unit 7,
Pocasset, MA 02559 USA

Phone: 508 563 1800

Email: geusa@geicp.com

Europe

Friedenbachstrasse 9,
35781 Weilburg,
Germany


Phone: +49 6471 3778517

Email: gegmbh@geicp.com

www.geicp.com




High Efficiency Sample Introduction System (HE-SIS)



HE-SIS Bracket Support

Every HE-SIS is designed to suit a specific instrument model, and includes an instrument-specific mounting bracket support.

Part Number	Description
KT-1155	HE-SIS for Agilent® ICP-MS
KT-1172	HE-SIS for TOFWERK icpTOF
KT-1172	HE-SIS for Thermo Scientific® Q, RQ, TQ ICP-MS
KT-1172	HE-SIS for Thermo Scientific® Neoma MC-ICP-MS
KT-1184	HE-SIS for PerkinElmer® NexION 1000, 1100, 2000, 2200, 5000 ICP-MS
KT-1204	HE-SIS for PerkinElmer® NexION 300, 350 ICP-MS
KT-1205	HE-SIS for NU ATTOM
KT-1213	HE-SIS for Thermo Scientific® X-Series
KT-1215	HE-SIS for Thermo Scientific® Neptune/Element
KT-1219	HE SIS for Nu Vitesse



HE-SIS for NexION 5000 ICP-MS

HE-SIS Kit Features



This specially designed concentric glass nebulizer is based on our popular MicroMist™ design, capable of efficiently nebulizing limited sample volumes at low sample and argon gas flow rates.



Our patent-pending MicroJet™ gas adapter shapes the nebulizer aerosol plume to reduce sample deposition on the spray chamber walls and enhance transport efficiency.



The Lotis™ HE spray chamber directly couples to the ICP-MS torch, providing the highest transport efficiency and excellent washout between samples.

Optimizing Operating Parameters




- Nebulizer gas flow rate (L/min)
- Nebulizer sample flow rate ($\mu\text{L}/\text{min}$)




- MicroJet gas flow rate (L/min)

*Combined gas flow rate through the injector is typically close to 1.0 L/min


Optimizing Nebulizer Gas Flow Rate



Optimizing MicroJet Gas Flow Rate



Average Sensitivity Ratio – Comparison Brand X



Demountable Torch – Interchangeable Injector




- A demountable torch provides the benefit of interchangeable injectors.
- 1.5mm and 2.5mm ID quartz studied
- Other injector ID's and materials available



* D-Torch™ for NexION 2200/5000 ICP-MS shown, P/N [30-808-3927](#)

Sensitivity Comparison – Injector ID



HE-SIS Literature

The HE-SIS has been coupled to many different ICP-MS platforms for a wide variety of applications, including single-cell, single-particle, nanoparticle, and low-volume sample studies, such as nanoplastics and microplastics with up to **95% transport efficiency**.

- *Development of single-cell ICP-TOFMS to measure nanoplastics association with human cells, Environ. Sci.: Nano, 2023, 10, 3439.*
- *Breaking barriers in Microplastic Detection using Single-Particle ICP-TOFMS, Lyndsey Hendriks, TOFWERK.*
- *Towards Automated Routine Analysis of the Distribution of Trace Elements in Single Cells using ICP-MS, Current Trends in Mass Spectrometry, March 2020.*
- *Very low mass isotope data collection with the Nu Vitesse, measurement of microplastic particles, Vitesse Note NT10.*
- In addition to many scientific presentations.

HE-SIS Summary

- In order to achieve optimum performance, it is necessary to optimize all operating conditions for both the instrument and sample introduction system.
- Our example showed the optimum sensitivity was observed at a nebulizer gas flow rate of 0.35 L/min and sample uptake in the range of 20 to 40 μ L/min.
- Glass Expansion's HE-SIS is 2–4x more sensitive than another commercially available system.
- Optimum make-up gas flow was dependent on the ID of the injector:
 - Smaller bore injector (1.5mm ID) provided highest sensitivity at a make-up gas flow of 0.50 L/min, combined gas flow of 0.85 L/min.
 - Larger bore injector (2.5mm ID) provided highest sensitivity at a make-up gas flow of 0.80 L/min, combined gas flow of 1.15 L/min.



Approaches and Strategies for the Detection and Quantification of Microplastics by Single Particle-ICP-MS



Contact Me

Aaron Hineman

**Atomic Spectroscopy Product Line
Leader, PerkinElmer**





Microplastics

Analytical Techniques

Microbeads in toiletries - Method 623.1

- Microbeads in products are extracted and analyzed to determine the composition of the beads. FT-IR is used to determine and provides confirmation of the presence or absence of plastic microbeads.



TGA-IR-GC/MS

- TGA - measures both the weight loss and the rate of evolution of products, which provides detailed kinetic information of polymers' decomposition.
- IR - quantitation of individual polymers and mixtures, functional group analysis of the volatile products.
- GC/MS - quantitation of polymers mixtures and functional group of the volatile products.




Microplastics

Analytical Techniques

Single Particle-ICP-MS

- Detection and counting of carbon-based particles
 - Short dwell times reduce the background
 - Accurate results for particles down to 2 μm
- Linear pass spray chamber is recommended to efficiently transports microplastics to the plasma
- SP-ICP-MS serves as a good screening technique for microplastics
 - Other techniques are required to determine composition



Carbon by Single Particle ICP-MS

Approaches and Strategies for Method Development

Platform selection

- ICP-MS vs ICP-MQ

Gas selection


- Std vs NH_3 vs. H_2

Isotope selection

- C_{12} vs C_{13}

Universal Cell Parameters

- Cell gas flow rate => Affect reaction rate
- Rejection parameter q (RPq) => Control product formation
- Rejection parameter a (RPa) => Attenuate selected element signal




Particle size
Detection
Limit

Instrumentation



NexION 5000



HE-SIS (P/N KT-1184)

Carbon by Single Particle ICP-MS

Operating Conditions

Parameter	Value
Neb Gas Flow (L/min)	0.46
AMS/Carrier Gas Flow (L/min)	0.6
Dwell Time (μ s)	25
IGM	Extraction
Scan Time (s)	30 - 180
Sample Flow Rate (μ L/min)	13
Transport Efficiency (%)	89.23

Carbon by Single Particle ICP-MS

Operating Conditions

Mode	Mass (amu)	Mode	Gas	Gas Flow	RPa	RPq	AFT	LOD (nm)
Standard	C13	MS/MS	-	-	0.017	0.25	150	1819
Standard	C12	MS/MS	-	-	0.02	0.25	150	990
DRC H ₂	C13	MS/MS	H ₂	1	-	0.45	225	672
DRC H ₂	C12	MS/MS	H ₂	1	0.043	0.45	200	722
DRC NH ₃	C13	MS/MS	NH ₃	0.1		0.45	125	1559
DRC NH ₃	C12	MS/MS	NH ₃	0.1	0.045	0.45	125	824

Carbon by Single Particle ICP-MS

Experimental

Standards

- 100 nm Au NP - N8151036
- Carbon 1000 ppm – Inorganic ventures

Samples

- Negative control - UPW
- 8 um polystyrene beads – positive control
- Tea bag – Sample
- Tap water – Sample blank



Carbon by Single Particle ICP-MS

Method Validation

Mode	Isotope	4 µm		8 µm	
		Size (nm)	Particle/mL	Size (nm)	Particle/mL
CRM Certified Values		4043	130000	7989	40000
Standard	C13	3426.8	226420*	9000.3	32579
Standard	C12	4035.3	123494	8723.1	31025
DRC H ₂	C13	3726.4	130561	8761.0	36371
DRC H ₂	C12	3515.3	119638	7636.9	35829
DRC NH ₃	C13	3135.6	123580	7834.5	44296
DRC NH ₃	C12	4125.9	133319	9095.2	42637

** Inaccurate counting due to low signal / high background*

Carbon by Single Particle ICP-MS


Results – 4 μm - C12 MS/MS DRC H_2



Average = 3.5 μm

Carbon by Single Particle ICP-MS


Results – 8 μm - C12 MS/MS DRC H_2



Average = 7.6 μm

Carbon by Single Particle ICP-MS

Mixed Particle Standards: 4 + 8 μm - C12 MS/MS DRC H₂



Carbon by Single Particle ICP-MS


Mixed Particle Standards: 4 + 8 μm - C12 MS/MS DRC H_2

Mode	Isotope	4 μm		8 μm	
		Size (nm)	Particle/mL	Size (nm)	Particle/mL
CRM Certified Values		4043	56500	7989	22500
Standard	C13	4234.2	44210	9046.1	24647
Standard	C12	4134.2	56189	8769.4	22407
DRC H_2	C13	4080.2	75406	9171.8	24992
DRC H_2	C12	3852.3	63543	8262.3	24072
DRC NH_3	C13	3968.6	44899	8271.0	24820
DRC NH_3	C12	4304.2	68943	9601.6	26543


Carbon by Single Particle ICP-MS

Tea Bag Analysis

- Brand 1




- Brand 2




Carbon by Single Particle ICP-MS

Tea Bag Analysis

- Brand 1



- Brand 2



Carbon by Single Particle ICP-MS

Tea Bag Analysis

Sample	Mode	amu	Start (nm)	End (nm)	Most Freq. Size (nm)	Mean Size (nm)	Part. Conc. (parts/mL)
Matrix Blk	DRC H ₂	C12	722	-	-	-	-
Tea Bag 1	DRC H ₂	C12	921	8671	3021	3568	3081
Tea Bag 2	DRC H ₂	C12	796	11146	3896	4549	10212



Carbon by Single Particle ICP-MS

Conclusions

- Carbon quantification by ICP-MS is feasible, various reaction gases, cell conditions, can be used to bring down the background to a measurable level.
 - Validated using DRC technology using 4 um and 8um standards to sub-micron detection limits
- Single particle can be used to screen for microplastics in a variety of matrices
- Improvement in detection limits is needed to achieve nanoplastics analysis
- Improvement in detection limits can be achieved with carbon free reagents and labware
- Microplastics counting can be achieved by SP-ICP-MS, but sizing is a challenge with unknown density
- Sample analysis time 2-3 min SP-ICP-MS vs 40-45 min using imaging FTIR.

More
Information
Online

Microplastics Analysis Solutions Spotlight


Resources

Applications

Blogs

Webinars

Technologies



Unlocking Carbon-13 with Single Particle ICP-MS: Feasibility Study for Microplastic Detection

Carbon is difficult to measure with ICP-MS because of its high ionization potential and its presence in both the argon used to generate the plasma (primarily in ...

[Learn More](#)



Customer Spotlight: Norwegian Institute for Water Research (NIVA)

Representatives from the Norwegian Institute for Water Research (NIVA) discuss their investigation into the presence of microplastics in the Arctic Ocean. The r ...

[Learn More](#)



Microplastics Analysis Brochure

PerkinElmer offers comprehensive and innovative solutions for the identification and quantification of microplastics. Our portfolio includes leading technologie ...

[Learn More](#)

Thank you, questions?

Aaron Hineman

Atomic Spectroscopy Product Line Leader



© 2025 PerkinElmer US LLC