

**Glass DC Nebulizers**

<b>Conikal</b>	
1mL/min uptake	A41-1-UC1S
<b>SeaSpray</b>	
2mL/min uptake	A41-1-USS2S
<b>Slurry</b>	
4mL/min uptake	A41-1-US6
<b>MicroMist</b>	
600µL/min uptake	A41-07-UM06S
200µL/min uptake	A41-07-UM02
100µL/min uptake	A41-07-UM01
50µL/min uptake	A41-07-UM005



**DC Nebulizers for HF**

<b>OpalMist</b>	
1mL/min uptake	A41-07-PFA1S
600µL/min uptake	A41-07-PFA06S
200µL/min uptake	A41-07-PFA02
100µL/min uptake	A41-07-PFA01
50µL/min uptake	A41-07-PFA005



<b>DuraMist</b>	
1mL/min uptake	A41-07-DM1S
400µL/min uptake	A41-07-DM04



<b>VeeSpray</b>	
0.6 - 3mL/min uptake	A41-1-CV6



**IsoMist XR**

<b>Programmable Temperature Cyclonic Spray Chamber</b>	
Glass spray chamber	KT-1067-XR
Quartz spray chamber	KT-1067Q-XR
PFA spray chamber	KT-1067P-XR



**Spray Chambers**

<b>Tracey</b>	
50mL Cyclonic	20-808-8881HE
With B14 gravity drain	



<b>Twister</b>	
50mL Cyclonic	20-809-9119HE
With B14 gravity drain	



<b>Twister</b>	
50mL Cyclonic	20-809-0222HE



<b>Cinnabar</b>	
20mL Cyclonic	20-809-0164HE



<b>Drain Device</b>	21-808-8214
---------------------	-------------



<b>Transfer Tube, Spray Chamber to ICPE-9000/9800 Torch</b>	
Glass Transfer Tube 86mm	21-809-2846



**Autosampler Probes** (See page 12 for complete list of probes)

**For Shimadzu ASC-6100F/9800**  
 PTFE Sheathed Carbon Fibre Probe  
 0.3mm ID with UniFit 70-803-1200

**For Shimadzu AS10**  
 PTFE Coated Vitcon O-rings  
 (PKT 10) With UniFit 70-803-1477

## Spray Chambers for HF

### Tracey TFE

50mL Cyclonic 20-809-2508  
 With B14 gravity drain  
 Suitable for use with HF



### HF Drain Device

21-809-2540



### Mounting bracket for HF Spray Chamber

21-809-2839



### Transfer Tube, Spray Chamber to ICPE-9000/9800 Torch

PFA Transfer Tube 82.5mm 31-808-3071

## Torches

### Semi Demountable Torch (HF)

30-808-8039

Suitable for ICPS-7500/8100 and ICPE-9000/9800

#### Comprises:

Torch Body 31-808-1258  
 Torch Adaptor 31-808-1260  
 Tapered Alumina Injector 1.5mm 31-808-0013  
 Joint Clip B14-Nickel KCM-14  
 GazFit Connectors (PKT 4) for 4mm OD side arms GAZ-04



### Tapered Injectors

Alumina Injector 1.0mm 31-808-0835  
 Alumina Injector 1.5mm 31-808-0013  
 Alumina Injector 1.8mm 31-808-0015  
 Alumina Injector 2.0mm 31-808-0027  
 Alumina Injector 2.4mm 31-808-0385



## Torches

### Capillary Injectors

Sapphire Injector 2.0mm 31-808-3001

### Quartz Torches


Quartz Torch Standard 1.5mm Injector (Suitable for 7500/8100) 30-808-1179

Quartz Torch for Organics, 1.0mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0519

Quartz Torch for High TDS, 1.8mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0518

Quartz Torch for High TDS, 2.0mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0564

Quartz Mini Torch 1.5mm Injector (Suitable for ICPE-9000/9800) 30-807-0534



## RF Coils

Copper 70-900-6000C  
 Silver 70-900-6000S  
 Gold 70-900-6000G

RF Coil Installation Kit 70-900-6000T

## Common Accessories

### Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1312



### Gaurdian In-Line Filters

Guardian in-line non-return gas filter 70-803-1362  
 Guardian in-line sample filter 70-803-1108



### Trident Internal Standard Kit

Standard Kit 60-808-1185  
 Standard Kit for HF 60-808-1150

