

## Glass DC Nebulizers

### SeaSpray

2mL/min uptake      A21-07-USS2

### MicroMist

600µL/min uptake      A21-07-UM06

### Slurry

4mL/min uptake      A21-07-US6

### Conikal

2mL/min uptake      A21-07-UC2



## Quartz DC Nebulizers

### SeaSpray

1mL/min uptake      A21-07-QSS1



## DC Nebulizers for HF

### DuraMist

1mL/min uptake      A21-07-DM1



### OpalMist

2mL/min uptake      A21-07-PFA2



## IsoMist XR

### Programmable Temperature Cyclonic Spray Chamber

IsoMist XR Kit      KT-1012-XR

#### Comprises:

Glass Twister Spray Chamber      21-809-3185

Assembly (LOL)      20-809-5175

IsoMist Mount      21-809-2773

Adaptor GAZ08-12 90deg      21-809-2871

Clamping Screw      21-808-0953

UniFit Right Angle Connector (PKT 3)      GAZ-04-90



## Spray Chambers

### Tracey

50mL Cyclonic      20-809-0434HE



### Twister

50mL Cyclonic (Aqueous)      20-809-0208HE  
With 7mm baffle



### Twister

50mL Cyclonic (Organics)      20-809-0446HE  
With 4mm baffle



### Twister

50mL Cyclonic (Aux. Port)      20-809-0004HE



### Cinnabar

20mL Cyclonic      20-809-9121HE



### HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

HydraMist Kit      KT-1162



## Spray Chambers for HF

### Tracey TFE

50mL Cyclonic      20-809-2510

Plugs directly into cassette mount



## Torches

**D-Torch** 30-808-3300 

**Comprises:**

- Base and Inner Tube 31-808-3292
- Quartz Outer Tube 31-808-3294
- Retaining Ring 31-808-2818

**D-Torch Options**

- Ceramic Outer Tube 31-808-2990 
- Copper Foils (PKT 12) for Cermic Outer Tube 31-808-3356
- O-Ring Kit 70-803-0762


**Quartz Torch**

- Without slots 31-808-3254 
- With single slot 31-808-3255
- With 3 slots 31-808-3256

**Quartz Torch Options:**

- Copper Foil for Quartz Torch or D-Torch Quartz Outer Tube (pack of 12) 31-808-3289

**Tapered Quartz Injectors**

- Quartz Injector 0.8mm 31-808-3257 
- Quartz Injector 2.0mm 31-808-3258
- Alumina Injector 1.2mm 31-808-3287
- Alumina Injector 2.0mm 31-808-3286 

**Sapphire Injectors**

- Sapphire Injector 2mm 31-808-3528

**Purge Windows**

- Axial 17mm disk 31-808-0551 
- Axial 32mm disk 31-808-2649
- Radial (Post Nov. 2004) 31-808-2788

## Common Accessories

**Elegra Argon Humidifier**

- Elegra Argon Humidifier 70-803-1269 

## Common Accessories

(continued)

**Guardian In-Line Filters**

- Guardian DC in-line non-return gas filter 70-803-1942 
- Guardian in-line sample filter 70-803-1108
- Eluo Adaptor for Inline Filter 70-803-1160

**Trident Internal Standard Kits**

- Trident CT Standard Kit 60-703-1179 
- Standard Kit for HF solutions 60-808-1150 

**Magnifier Inspection Tool** 70-803-1923 

## Autosampler Probes (See page 22-23 for complete list of probes)

**Suitable for PE S10 or AS93+**

- Guardian Probe 0.75mm Probe Connecting Line (Red) (\*for PE S10 only) 70-803-2851 
- Guardian Probe 0.75mm Probe Connecting Line (Red) (\*for AS93+ only) 70-803-2836 
- PTFE Encapsulated Carbon Fibre Probe 0.25mm ID 70-803-1071 
- PTFE Encapsulated Carbon Fibre Probe 0.5mm ID 70-803-1440
- PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0816

**Suitable for Cetac ASX-200/500/800 Series**

- Guardian Probe 0.75mm Probe Connecting Line (Red) 70-803-1803 
- Polyimide sheathed Autosampler Probe 0.5mm ID 60-808-1186L
- PTFE Encapsulated Carbon Fibre Probe 0.5mm ID 70-803-0784 
- PTFE Encapsulated Carbon Fibre Probe 0.75mm ID 70-803-1443
- PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0793

**For Cetac ASX-7400, 7600 Series**

- Guardian Probe 0.75mm Probe Connecting Line (Red) 70-803-2940 