
## Eluo Nebulizer Cleaning Tool

As easy as 1, 2, 3



Particle build-up in a nebulizer capillary and tip causes sample flow to be constricted, thereby reducing nebulizer efficiency and performance. Now, blocked nebulizers can be safely and easily restored to optimum performance with the Eluo.

The Eluo is designed to efficiently deliver a cleaning solution through the nebulizer capillary to dislodge particle build-up and thoroughly clean the nebulizer. One simple action does it all. No more messy procedures or shattered nebulizers in ultrasonic baths. Use the Eluo regularly to maintain nebulizer performance and prolong nebulizer life. Every lab should have an Eluo.

Order P/N <u>70-ELUO</u> for glass concentric nebulizers or P/N <u>70-ELUO-OPD</u> for OpalMist or DuraMist nebulizers.

Ordering Information	
Part Number	Description
70-ELUO	Eluo Nebulizer Cleaning Tool (For glass concentric nebulizers)
70-ELUO-OPD	Eluo Nebulizer Cleaning Tool (For OpalMist or DuraMist nebulizers)



P/N 70-ELUO-OPD

For more information visit: www.geicp.com



## Asia Pacific

6 Central Boulevard, Port Melbourne, Vic 3207, Australia

(61) 3 9320 1111 enquiries@geicp.com

## **Americas**

31 Jonathan Bourne Drive, Unit 7 Pocasset, MA 02559, USA

508 563 1800 geusa@geicp.com

## **Europe**

Friedenbachstrasse 9, 35781 Weilburg, Germany +49 6471 3778517

gegmbh@geicp.com